

Combining analysis of solid and liquid biopsy: insights into tumor heterogeneity

by **L Alvarez^{1,2,3}, A Rodrigo^{1,2}, A Terradez^{1,2}, A Finzel Pérez¹, JF Laes¹**

1: OncoDNA SA, Gosselies, Belgium

2: Biosequence SL, Valencia, Spain

3: I3S – Instituto de Investigação e Inovação da Universidade do Porto, Porto, Portugal

l.alvarez@oncodna.com, a.rodrigo@oncodna.com, a.terradez@oncodna.com, a.finzel-perez@oncodna.com, jf.laes@oncodna.com

INTRODUCTION

Liquid biopsies are a promising tool for molecular tumor profiling and monitoring (Siravegna *et al.*, 2017). One of their main advantages over solid biopsies is that they can detect intratumoral heterogeneity, which is unappreciated by single-site tissue biopsies. However, we are still at the beginning of their incorporation into routine oncology practice.

Aim of the project: *What is the clinical utility of combining solid and liquid biopsies in the molecular profiling of tumor patients?*

MATERIALS AND METHODS

We analysed 112 samples from patients with different cancer types, stage III or IV, using **OncoSTRAT&GO™** (OncoDNA):


- **Solid biopsy (FFPE block)** 
 - Next-generation sequencing → more than 200 genes
 - (Immunohistochemistry)
- **Liquid biopsy (blood sample)** 
 - Next-generation sequencing of ctDNA → **27 genes**

RESULTS


First, we wanted to study the concordance between solid and liquid biopsies in detecting actionable variants present in the tumors → **60.7% complete concordance** (68 out of 112 samples):

TUMOR TYPE	PATIENT NUMBER (%)	CONCORDANCE	
		PATIENT NUMBER (%)	
Breast Cancer	29 (25.9)	14 (48.3)	
Lung Cancer	15 (13.4)	11 (73.3)	
Colorectal Cancer	13 (11.6)	8 (61.5)	
Gynecologic Cancer	12 (10.7)	5 (41.7)	
Pancreatic Cancer	9 (8.0)	5 (55.6)	
Others	34 (30.4)	25 (73.5)	
TOTAL	112	68 (60.7)	

The rest of the samples (39.3%) presented different discrepant patterns (in orange tonalities):



We then compared the concordant and discrepant variant allele frequencies, and they showed similar distributions with no significant statistical differences: mean values of 7.8/9.4% (P= 0.58, Mann-Whitney test) and 34.9/25.7% (P= 0.08, Mann-Whitney test) in liquid and solid biopsy, respectively:



CASE STUDY

Colorectal cancer patient

 **TP53 (L194R, 39.34%)** → No related treatments available

 **TP53 (L194R, 33.78%)**
KRAS (G12D, 13.13%)
PIK3CA (E545A, 12.24%) → **KRAS (G12D)**: Metastatic patients should not be treated with cetuximab or panitumumab (FDA)

PIK3CA (E545A): PI3K/mTOR dual inhibitor might be useful (ASCO 2015)

CONCLUSIONS

The discrepancy found in the data cannot be explained by the sensibility of the analysis, fact that points to a key influence of the heterogeneity of the tumor, which is only captured by liquid biopsies. So, combining solid and liquid biopsies would give a full picture of the tumor.

Indeed, our results show the usefulness of the combination of both biopsies in clinical practice, as it provides additional information in 39.3% of the cases. This leads to a broader characterization of the tumor molecular profile and therefore to a better disease management.